













# CAPSULE POLISHING

Designed for polishing and cleaning the residual cortex from the capsular bag. Silicone tipped cannulae provide gentle polishing and atraumatic insertion and removal. Micro etched cannulae are slightly roughened for gentle polishing.


## Micro Etched

Code	Gauge and size	Description	Tip
M8604	 21G x 1" 0.81 x 25mm	Capsule polisher. Angled 12.5mm from micro etched tip. Aspiration through 0.3mm side port	
M8604A	 23G x 7/8" 0.64 x 22mm	Capsule polisher. Angled 9mm from micro etched tip. Aspiration through 0.3mm side port	
M8604B	 21G x 7/8" 0.81 x 22mm	Capsule polisher. Angled 9mm from micro etched tip. Aspiration through 0.3mm side port	
M8603	 21G x 7/8" 0.81 x 22mm	Capsule polisher. Angled 3mm from end. Micro etched 5mm from end	


## Silicone Tipped

Code	Gauge and size	Description	Tip
M8605	 27G x 7/8" 0.4 x 22mm	Silicone tipped polisher	
M8605A	 23G x 7/8" 0.64 x 22mm	Silicone tipped polisher	
M8608	 27G x 7/8" 0.4 x 22mm	Flexible tipped polisher. Curved through 45°	

## Posterior Scraper

Code	Gauge and size	Description	Tip
M8645	 23G x 7/8" 0.64 x 22mm	Posterior Capsule scraper	

## Olive Tipped

Code	Gauge and size	Description	Tip
M8620	 23G x 1" 0.64 x 25mm	Olive Tipped Polisher. Straight	
M8620C	 23G x 1" 0.64 x 25mm	Olive Tipped Polisher. Curved	
M8621	 23G x 1" 0.64 x 25mm	Olive Tipped Polisher. Curved, micro etched tip	
M8622	 25G x 7/8" 0.5 x 22mm	Olive Tipped Polisher. Angled, micro etched tip	