




STERIMEDIX[®]

single-use ophthalmic innovations





delivering **innovation** in ophthalmic surgery



continuous innovation in response to our customers' requirements; continued growth from excellent product quality and customer service.

This catalogue contains details of the entire Sterimedix range, including unique new products designed to work with the latest techniques, in addition to a tried and trusted portfolio of single use cannulae for all ophthalmic procedures.



Sterimedix. Dedicated to supporting our distributors and their customers

Our quest is to continually improve both our service and products, supplying innovative single use ophthalmic instruments at competitive prices, presented in user friendly packaging, and meeting our customers' delivery requirements.

Sterimedix products are packed sterile in a single blister in space saving boxes designed to minimise waste. Our products may also be supplied non-sterile in formats compatible with the requirements of kit packers.

Single use products

Page	4 - 5	NEW from Sterimedix	15	Irrigating / Aspirating Cannulae
	6	Refractive / Scleral Marker	16 - 17	Bimanual Handles and Cannulae
	7	Anaesthesia	18	I/A Handpieces
	8	Anterior Chamber Maintenance	19	Capsule Polishing
	9	Viscoexpression / Viscoelastic	20	Lachrymal Cannulae
	10 - 11	Capsulotomy	21 - 23	Vitreoretinal Products
	12 - 13	Hydrodissection / Hydrodelineation	24	Drapes and Fluid Collection
	14	Nucleus Removal / Cortex Removal	25 - 27	Product Index
			27	Product Key / Kit Packers

QUALITY ASSURANCE Sterimedix Limited is committed to producing products of the highest quality. The company's quality systems and practices are regularly audited and have ISO 13485:2003, ISO 9001:2000 and CE accreditations together with United States FDA registration, and the appropriate approvals in all other markets where Sterimedix is active.

NEW FROM STERIMEDIX

Continuing our reputation for innovation, we have introduced a series of new cannulae to complement the existing range. In answer to the needs of our customers, we continue to develop new products to match changing techniques and surgical advances.

REFRACTIVE



- NEW M2823 ● 25G x 7/8" 0.5 x 22mm Irrigating Lasik/Intralasik Cannula (Daya).
Formed 10.5mm from flattened spatulated tip
- NEW M2824 ● 23G x 7/8" 0.64 x 22mm Irrigating Lasik/Intralasik Cannula (Daya).
Formed 10.5mm from flattened spatulated tip



ANAESTHESIA



- NEW M0279 ● 20G x 1" 0.9 x 25mm Sub Tenon Cannula. Curved with Flattened Tip



ANTERIOR CHAMBER MAINTENANCE



- NEW M2270A ● 30G x 7/8" 0.3 x 22mm Anterior Chamber (Rycroft) Angled 35°, 8mm from tip
- NEW M2279 ● 27G x 7/8" 0.4 x 22mm Anterior Chamber Angled 45°, 11mm from tip



VISCOELASTIC



- NEW M3917 ● 27G x 7/8" 0.4 x 22mm Viscoelastic. Formed 11mm from tip



- NEW M3927 ● 27G x 7/8" 0.4 x 22mm Viscoelastic. Angled at 25°, 8mm from tip

HYDRODISSECTION



NEW M2273F



25G x 7/8"

0.5 x 22mm

Helsinki Hydrodissection Cannula



NEW M3912



23G x 7/8"

0.64 x 22mm

Hydrodissector (Mackool) Angled 30°, 8mm from tip and flattened 4mm from tip



NEW M3913



25G x 7/8"

0.5 x 22mm

Hydrodissector (Daya), Formed 10.5mm from flattened tip



NEW M3914



23G x 7/8"

0.64 x 22mm

Hydrodissector (Daya), Formed 10.5mm from flattened tip

BIMANUAL HANDLES



NEW M5513




23G x 5/8"

0.64 x 16mm

Cannula Thin Wall. Curved with twin 0.4mm side ports



BIMANUAL I/A SYSTEMS



NEW M5535



23G x 5/8"

0.64 x 16mm

Irrigating Cannula Thin Wall. Curved with twin 0.4mm side ports



23G x 5/8"

0.64 x 16mm

Aspirating Cannula Thin Wall. Curved with 0.3mm side port and micro etched tip



VITREORETINAL



NEW M4470C



20G/23G 1 1/4"

















0.9/0.64 x 32mm

Curved Reducing Cannula. Reduced 3mm from end of tip



REFRACTIVE

Products designed for laser refractive procedures.

Code	Gauge and size	Description	Tip
M2804	 25G x 7/8" 0.5 x 22mm	Irrigating Lasik Cannula. Formed 9mm from the tip. Six 0.15mm ports at 2, 4 and 6mm from tip	
M2817	 27G x 7/8" 0.4 x 22mm	Irrigating Lasik Cannula. Formed 9mm from the tip. Four 0.15mm ports at 2 and 4mm from tip	
M2819	 27G x 7/8" 0.4 x 22mm	Irrigating Lasik Cannula. Angled at 25°, 7mm from the tip.	
M2822	 25G x 7/8" 0.5 x 22mm	Irrigating Lasik Cannula. Formed 9mm from the tip. Six 0.15mm ports at 2, 4 and 6mm from closed tip	
M2801	 25G x 7/8" 0.5 x 22mm	Irrigating Lasik Cannula (Lin). Angled with two 0.15mm ports and flattened tip	
M2802	 25G x 1 1/8" 0.5 x 29mm	Irrigating Lasik Cannula (Buratto). Formed with two 0.15mm side ports	
NEW M2823	 25G x 7/8" 0.5 x 22mm	Irrigating Lasik/Intralasik Cannula (Daya). Formed 10.5mm from flattened spatulated tip	
NEW M2824	 23G x 7/8" 0.64 x 22mm	Irrigating Lasik/Intralasik Cannula (Daya). Formed 10.5mm from flattened spatulated tip	

SCLERAL MARKER

The scleral marker is designed to be used as part of the intravitreal injection procedure for the treatment of AMD. It incorporates two calibrated markers to mark the injection point.




Code	Size	Description
M7500	3.5mm / 4.0mm	Scleral Marker. 80mm long





ANAESTHESIA

Retrobulbar

For injecting anaesthetic into the muscle cone.



Code	Gauge and size	Description	Tip	
 M0275	● 25G x 1½"	0.5 x 38mm	Retrobulbar Tri-Facet	
M0637	● 25G x 1½"	0.5 x 38mm	Retrobulbar (Atkinson)	
M0638	● 23G x 1½"	0.64 x 38mm	Retrobulbar Thin Wall (Atkinson)	

Facial Infiltration

Code	Gauge and size	Description	Tip	
 M0285	● 30G x 5/8"	0.3 x 16mm	Facial Infiltration	



Sub Tenon Cannula

For injecting anaesthetic into the sub tenon space.

Code	Gauge and size	Description	Tip	
 M0278	● 19G x 1"	1.1 x 25mm	Sub Tenon Cannula. Curved with Flattened Tip	
NEW M0279	● 20G x 1"	0.9 x 25mm	Sub Tenon Cannula. Curved with Flattened Tip	

Peribulbar

For injecting anaesthetic around the globe.

Code	Gauge and size	Description	Tip	
 M0641A	● 27G x 1¼"	0.4 x 32mm	Peribulbar	
M0642	● 25G x 1¼"	0.5 x 32mm	Peribulbar	
M0643	● 23G x 1¼"	0.64 x 32mm	Peribulbar	
M0641	● 27G x 7/8"	0.4 x 22mm	Peribulbar	
M0642A	● 25G x 7/8"	0.5 x 22mm	Peribulbar	







ANTERIOR CHAMBER MAINTENANCE

For maintaining and forming the anterior chamber during surgery by injecting or removing viscoelastic materials, BSS or air.






Anterior Chamber Cannulae

Code	Gauge and size	Description	Tip
			
M2639	● 25G x 1" 0.5 x 25mm	Anterior Chamber (Rycroft) Angled 35°, 7.5mm from tip	
M2270	● 30G x 7/8" 0.3 x 22mm	Anterior Chamber (Rycroft) Angled 35°, 4mm from tip	
NEW M2270A	● 30G x 7/8" 0.3 x 22mm	Anterior Chamber (Rycroft) Angled 35°, 8mm from tip	
M2273	● 30G x 7/8" 0.3 x 22mm	Anterior Chamber (Rycroft) Angled 35°, 5.5mm from tip	
M2273A	● 27G x 7/8" 0.4 x 22mm	Anterior Chamber (Rycroft) Angled 35°, 4.5mm from tip	
NEW M2279	● 27G x 7/8" 0.4 x 22mm	Anterior Chamber Angled 45°, 11mm from tip	
M2275	● 25G x 7/8" 0.5 x 22mm	Anterior Chamber (Rycroft) Angled 35°, 5mm from tip	
M2273G	● 23G x 7/8" 0.64 x 22mm	Anterior Chamber (Rycroft) Angled 35°, 5.5mm from tip	
M2642	● 21G x 7/8" 0.81 x 22mm	Anterior Chamber (Rycroft) Angled 35°, 7.5mm from tip	
M2273B	● 20G x 7/8" 0.9 x 22mm	Anterior Chamber (Rycroft) Angled 45°, 7.5mm from tip	

Irrigating Cannulae




Code	Gauge and size	Description	Tip
			
M9280	● 19G x 1" 1.1 x 25mm	Washout Cannula (Bishop Harmon) Angled 32°, 8mm from tip	
			
M9281	○ 16G x 1" 1.6 x 25mm	Washout Cannula. Angled 32°, 6.5mm from tip	

Anterior Chamber Maintainers







Code	Gauge and size	Description	Tip
			
M2656	● 20G x 9/64" 0.9 x 3.5mm	Anterior Chamber Maintainer (Dahan) 150mm Silicone Tubing Micro etched allowing cannula to be self-retaining	
			
M2658	● 23G x 9/64" 0.64 x 3.5mm	Anterior Chamber Maintainer (Lewicky) 150mm Silicone Tubing Annular rings allow cannula to be self-retaining	
			
M2657	● 20G x 9/64" 0.9 x 3.5mm	Anterior Chamber Maintainer (Lewicky) 150mm Silicone Tubing Annular rings allow cannula to be self-retaining	

VISCOEXPRESSION/VISCOELASTIC

Viscoexpression Cannulae

Code	Gauge and size	Description	Tip
			
<p>The shape allows the cannulae to be inserted under the nucleus for expression using viscoelastic.</p>			
M3899	● 25G x 7/8" 0.5 x 22mm	Corydon. Flat tip angled 2.5mm from end	
M3899A	● 25G x 7/8" 0.5 x 22mm	Corydon. Rounded tip angled 1.8mm from end	

Viscoelastic Cannulae

Code	Gauge and size	Description	Tip
			
<p>This range of cannulae has smooth rounded tips for the safe introduction of viscoelastic.</p>			
M3926A	● 27G x 7/8" 0.4 x 22mm	Viscoelastic. Angled at 35°, 9mm from tip	
NEW M3927	● 27G x 7/8" 0.4 x 22mm	Viscoelastic. Angled at 25°, 8mm from tip	
M2272	● 26G x 7/8" 0.45 x 22mm	Viscoelastic. Angled at 35°, 9mm from tip	
M2639A	● 25G x 7/8" 0.5 x 22mm	Viscoelastic. Angled at 35°, 9mm from tip	
M3909	● 23G x 7/8" 0.64 x 22mm	Viscoelastic. Angled at 45°, 9mm from tip	
M3910A	● 21G x 7/8" 0.81 x 22mm	Viscoelastic. Angled at 45°, 9mm from tip	
			
M3910	● 21G x 7/8" 0.81 x 22mm	Viscoelastic (Straight)	
			
NEW M3917	● 27G x 7/8" 0.4 x 22mm	Viscoelastic. Formed 11mm from tip	



















CAPSULOTOMY

Irrigating cystotomes are designed for capsularhexis, intercapsular/endocapsular and can opener capsulotomy techniques.

Formed cystotomes are shaped to conform to the crystalline lens. Short curved cystotomes are more suited to deep set eyes or small pupils.







Capsularhexis cystotomes have a special tip which allows greater control of the flap when performing continuous curvilinear capsulotomies.

Irrigating Cystotomes



Code		Gauge and size		Description	Tip
	M3630	 30G x 5/8"	0.3 x 16mm	Cystotome (Formed)	
	M3610B	 27G x 5/8"	0.4 x 16mm	Cystotome (Formed)	
	M3610	 25G x 5/8"	0.5 x 16mm	Cystotome (Formed)	
	M3610S	 25G x 5/8"	0.5 x 16mm	Cystotome (Formed) Short Tip	
	M3602	 25G x 5/8"	0.5 x 16mm	Side Cutting Cystotome (Formed)	
	M3610C	 25G x 5/8"	0.5 x 16mm	Cystotome (Formed) Tip angled 45°	
	M3627A	 27G x 5/8"	0.4 x 16mm	Cystotome (Formed) Short Curve	
	M3627	 25G x 5/8"	0.5 x 16mm	Cystotome (Formed) Short Curve	
	M3615	 27G x 7/8"	0.4 x 22mm	Long Cystotome (Formed)	
	M3615A	 25G x 7/8"	0.5 x 22mm	Long Cystotome (Formed)	






Capsularhexis Cystotomes

Code	Gauge and size	Description	Tip
 M3000F	● 27G x 5/8" 0.4 x 16mm	Cystotome (Formed)	
M3000G	● 25G x 5/8" 0.5 x 16mm	Cystotome (Formed)	
 M3000E	● 27G x 5/8" 0.4 x 16mm	Cystotome (Formed) Short Curve	
M3000D	● 25G x 5/8" 0.5 x 16mm	Cystotome (Formed) Short Curve	




Curved Cystotomes

Code	Gauge and size	Description	Tip
 M3001	● 27G x 5/8" 0.4 x 16mm	Endocapsular Cystotome (Joyce)	

Angled Cystotomes

Code	Gauge and size	Description	Tip
 M3626	● 27G x 5/8" 0.4 x 16mm	Cystotome	
M3626A	● 25G x 5/8" 0.5 x 16mm	Cystotome	

















Reverse Formed

Code	Gauge and size	Description	Tip
 M3629	● 25G x 5/8" 0.5 x 16mm	Cystotome (Formed, Reverse Tip)	
M3629A	● 27G x 5/8" 0.4 x 16mm	Cystotome (Formed, Reverse Tip)	







HYDRODISSECTION

A range of cannulae with tips designed to deliver fluids to facilitate the separation of the cortex from the nucleus and capsule.





Hydrodissection Cannulae

Code	Gauge and size	Description	Tip
		Tip flattened in vertical plane.	
M2273C	● 27G x 7/8" 0.4 x 22mm	Sauter Hydrodissector (Curved)	
M2273D	● 25G x 7/8" 0.5 x 22mm	Sauter Hydrodissector (Curved)	
		Tip flattened in horizontal plane.	
M2273E	● 27G x 7/8" 0.4 x 22mm	Helsinki Hydrodissector (Curved)	
NEW M2273F	● 25G x 7/8" 0.5 x 22mm	Helsinki Hydrodissection Cannula	
		Tip flattened in vertical plane.	
M2274	● 25G x 7/8" 0.5 x 22mm	Kellan Hydrodissector (Curved)	
		Tip flattened in horizontal plane.	
M3992	● 30G x 7/8" 0.3 x 22mm	Hydrodissector (Angled 7mm from tip)	
M3902	● 27G x 7/8" 0.4 x 22mm	Hydrodissector (Angled 9mm from tip)	
M3900	● 25G x 7/8" 0.5 x 22mm	Hydrodissector (Angled 9mm from tip)	
			
M2273H	● 27G x 7/8" 0.4 x 22mm	Hydrodissector (Curved) with rounded tip	
			
M3904	● 27G x 3/4" 0.4 x 18mm	Hydrodissector (Angled 90°) (Chang)	
		Tip flattened in horizontal plane.	
NEW M3912	● 23G x 7/8" 0.64 x 22mm	Hydrodissector (Mackool) Angled 30°, 8mm from tip and flattened 4mm from tip	
		Tip flattened in horizontal plane.	
NEW M3913	● 25G x 7/8" 0.5 x 22mm	Hydrodissector (Daya), Formed 10.5mm from flattened tip	
NEW M3914	● 23G x 7/8" 0.64 x 22mm	Hydrodissector (Daya), Formed 10.5mm from flattened tip	



Cortical Cleaving Cannulae

Code	Gauge and size	Description	Tip
M3908	 25G x 7/8" 0.5 x 22mm	Cortical Cleaving Hydrodissector. Tapered to 0.30mm at tip and angled at 45°, 1mm from tip	
M3911	 27G x 7/8" 0.4 x 22mm	Cortical Cleaving Hydrodissector. Angled 11mm from tip with tip flattened in horizontal plane 2.75mm from end	
M3911A	 27G x 7/8" 0.4 x 22mm	Cortical Cleaving Hydrodissector. Angled 11mm from round tip	

J-Shaped Hydrodissection Cannulae







Code	Gauge and size	Description	Tip
M3901	 25G x 7/8" 0.5 x 22mm	J-Shaped Hydrodissector. Tapered and flattened tip, 1.75mm wide	
M3906A	 25G x 7/8" 0.5 x 22mm	J-Shaped Hydrodissector. Round tip, 1.75mm wide	

Hydrodelineation Cannulae

Code	Gauge and size	Description	Tip
M3907	 25G x 7/8" 0.5 x 22mm	Hydrodelineator. Tapered to 0.30mm at tip and angled at 45°, 4mm from tip	

















NUCLEUS REMOVAL

Irrigating Vectis designed for irrigating and removing the nucleus and the larger pieces of cortical debris whilst maintaining the anterior chamber depth.

Code	Gauge and size	Description	Tip
M4619	 25G x 1½" 0.5 x 38mm	Irrigating Vectis. 5mm wide, single port	
M4620	 25G x 1½" 0.5 x 38mm	Irrigating Vectis. 5mm wide, twin ports	
M4621	 27G x 1½" 0.4 x 38mm	Small Incision Irrigating Vectis. 3mm x 2.75 wide, single port	


CORTEX REMOVAL

For aspirating cortical debris when chamber is maintained with an anterior chamber maintainer or viscoelastic substance. The 'J' shaped cannulae are designed to retract the iris and aspirate difficult to reach cortex between the 10 to 2 O'clock position.

Code	Gauge and size	Description	Tip
M6576	 23G x 7⁄8" 0.64 x 22mm	Cortex Extractor. 0.3mm top port	
M6575	 23G x 5⁄8" 0.64 x 16mm	Cortex Extractor. 0.3mm top port	
M6575A	 23G x 5⁄8" 0.64 x 16mm	Cortex Extractor. 0.4mm top port	
M6576A	 23G x 1⁄2" 0.64 x 12.5mm	Cortex Extractor. 0.3mm top port	
M6525	 21G x 5⁄8" 0.81 x 16mm	Cortex Extractor. 0.3mm top port	
M6594	 23G x 5⁄8" 0.64 x 16mm	Cortex Extractor (Charleux). 0.3mm bevelled end port	
M6596	 24G x 1⁄2" 0.55 x 12.5mm	Cortex Extractor (Charleux). 0.3mm bevelled end port	
M6595	 23G x 1⁄2" 0.64 x 12.5mm	Cortex Extractor (Charleux). 0.3mm bevelled end port	
M5620	 23G x 1" 0.64 x 25mm	Binkhorst Hook (Right). 2.5mm wide 'J' shape	
M5622	 25G x 7⁄8" 0.5 x 22mm	Binkhorst Hook (Right). 2.0mm wide 'J' shape	
M5621	 23G x 1" 0.64 x 25mm	Binkhorst Hook (Left). 2.5mm wide 'J' shape	
M5623	 25G x 7⁄8" 0.5 x 22mm	Binkhorst Hook (Left). 2.0mm wide 'J' shape	

IRRIGATING/ASPIRATING CANNULAE

For removing cortical debris or viscoelastic substances while maintaining the anterior chamber during extracapsular and intercapsular/endocapsular cataract procedures.



Code	Gauge and size	Description	Tip
M6600	● 23G/23G 0.64/0.64mm	I & A Cannula (Standard). 0.4mm port (Simcoe) 200mm silicone tubing 2.65 OD x 0.65 ID	
M6601	● 23G/23G 0.64/0.64mm	I & A Cannula (Reverse). 0.4mm port (Simcoe) 200mm silicone tubing 2.65 OD x 0.65 ID	
M6605	● 23G/23G 0.64/0.64mm	I & A Cannula Short Pattern (Standard). 0.3mm port (Simcoe) 200mm silicone tubing 2.00 OD x 1.00 ID	

I/A Cannulae

Used for irrigating/aspirating within the anterior chamber.

Code	Gauge and size	Description	Tip
M6611	● 21G x 1½" 0.81 x 38mm	Straight Cannula	

Coaxial I/A Cannulae




Code	Gauge and size	Description	Tip
M6634	18G/21G 1.24/0.81mm	Coaxial I & A Cannula. 0.3 side aspiration port on inner cannula	
M6634B	18G/21G 1.24/0.81mm	Coaxial I & A Cannula. 0.3 side aspiration port on inner cannula. Tip angled 45°	









BIMANUAL HANDLES

All Bimanual handles have a coloured indicator to show the gauge size in accordance with ISO 6009.

Aspiration

Code	Gauge and size	Description	Tip
			
M5501	● 23G x 5/8"	0.64 x 16mm Cannula. Curved with 0.4mm side port and micro etched tip	
M5502	● 23G x 5/8"	0.64 x 16mm Cannula. Curved with 0.3mm side port and micro etched tip	
M5512	● 21G x 5/8"	0.81 x 16mm Cannula. Curved with 0.3mm side port and micro etched tip	
M5520	● 21G x 5/8"	0.81 x 16mm Cannula Thin Wall. Curved with 0.4mm side port and micro etched tip	
M5511	● 20G x 3/4"	0.90 x 19mm Cannula. Curved with 0.3mm side port and micro etched tip	
M5518	● 23G x 5/8"	0.64 x 16mm Cannula. Curved with 0.3mm side port	
M5515	● 21G x 7/8"	0.81 x 22mm Cannula. Curved with 0.3mm side port	



















Irrigation

Code	Gauge and size	Description	Tip
			
M5519	● 23G x 5/8"	0.64 x 16mm Cannula Ultra Thin Wall. Curved with 45° bevelled tip	
M5505	● 21G x 5/8"	0.81 x 16mm Cannula. Curved with 45° bevelled tip	
M5510	● 20G x 3/4"	0.90 x 19mm Cannula. Curved with twin 0.5mm side ports	
M5514	● 21G x 5/8"	0.81 x 16mm Cannula. Curved with twin 0.5mm side ports	
M5521	● 21G x 5/8"	0.81 x 16mm Cannula Thin Wall. Curved with twin 0.6mm side ports	
NEW M5513	● 23G x 5/8"	0.64 x 16mm Cannula Thin Wall. Curved with twin 0.4mm side ports	
M5523	● 23G x 5/8"	0.64 x 16mm Cannula Ultra Thin Wall. Curved with open end	
M5522	● 21G x 5/8"	0.81 x 16mm Cannula Thin Wall. Curved with open end	
			
M5509	● 21G x 5/8"	0.81 x 16mm Cannula. Straight with twin 0.5mm side ports	






Bimanual I/A Systems

A combination set of one Irrigating Handle and one Aspirating Handle manufactured using thin wall stainless steel tubing.










Code	Gauge and size	Description	Tip	
M5531	 21G x 5/8"	0.81 x 16mm	 	 
	 21G x 5/8"	0.81 x 16mm		
NEW M5535	 23G x 5/8"	0.64 x 16mm	 	 
	 23G x 5/8"	0.64 x 16mm		
M5532	 21G x 5/8"	0.81 x 16mm	 	 
	 21G x 5/8"	0.81 x 16mm		

For more information of other systems available, please contact Sterimedix or your local distributor.








Bimanual Aspirating Cannulae

Code	Gauge and size	Description	Tip
M6524	 21G x 5/8"	0.81 x 16mm	 

Bimanual Irrigating Cannulae

Code	Gauge and size	Description	Tip
M5405	 21G x 5/8"	0.81 x 16mm	 
M5414	 21G x 5/8"	0.81 x 16mm	 
M5409	 21G x 5/8"	0.81 x 16mm	 

COAXIAL I/A HANDPIECES









Code	Gauge and size	Description	Tip
			
M5580	16/21G 1.6/0.81mm	I/A Handpiece. Straight tip	
M5582	16/21G 1.6/0.81mm	I/A Handpiece. 45° angled tip	
M5583	16/21G 1.6/0.81mm	I/A Handpiece. 60° angled tip	
M5584	16/21G 1.6/0.81mm	I/A Handpiece. 90° angled tip	
			
M5581	16/21G 1.6/0.81mm	I/A Handpiece. Curved tip	










CAPSULE POLISHING

Designed for polishing and cleaning the residual cortex from the capsular bag. Silicone tipped cannulae provide gentle polishing and atraumatic insertion and removal. Micro etched cannulae are slightly roughened for gentle polishing.

Micro Etched

Code	Gauge and size	Description	Tip
 M8604	 21G x 1" 0.81 x 25mm	Capsule polisher. Angled 12.5mm from micro etched tip. Aspiration through 0.3mm side port	
M8604A	 23G x 7/8" 0.64 x 22mm	Capsule polisher. Angled 9mm from micro etched tip. Aspiration through 0.3mm side port	
M8604B	 21G x 7/8" 0.81 x 22mm	Capsule polisher. Angled 9mm from micro etched tip. Aspiration through 0.3mm side port	
 M8603	 21G x 7/8" 0.81 x 22mm	Capsule polisher. Angled 3mm from end. Micro etched 5mm from end	

Silicone Tipped

Code	Gauge and size	Description	Tip
 M8605	 27G x 7/8" 0.4 x 22mm	Silicone tipped polisher	
M8605A	 23G x 7/8" 0.64 x 22mm	Silicone tipped polisher	
 M8608	 27G x 7/8" 0.4 x 22mm	Flexible tipped polisher. Curved through 45°	

Posterior Scraper

Code	Gauge and size	Description	Tip
 M8645	 23G x 7/8" 0.64 x 22mm	Posterior Capsule scraper	

Olive Tipped

Code	Gauge and size	Description	Tip
 M8620	 23G x 1" 0.64 x 25mm	Olive Tipped Polisher. Straight	
 M8620C	 23G x 1" 0.64 x 25mm	Olive Tipped Polisher. Curved	
M8621	 23G x 1" 0.64 x 25mm	Olive Tipped Polisher. Curved, micro etched tip	
 M8622	 25G x 7/8" 0.5 x 22mm	Olive Tipped Polisher. Angled, micro etched tip	

LACHRYMAL CANNULAE

Cannulae for probing and irrigating the lachrymal ducts

Code	Gauge and size	Description	Tip
M9275	26G x 3/4"	0.45 x 20mm Lachrymal Cannula. Short Straight	
M9276	26G x 1 1/4"	0.45 x 32mm Lachrymal Cannula. Curved	
M9277	25G x 1 1/8"	0.5 x 29mm Lachrymal Cannula. Curved	
M9274	21G x 1 1/2"	0.81 x 38mm Lachrymal Cannula. Curved	
M9272	23G x 1/2"	0.64 x 12.5mm Lachrymal Cannula. Short Straight. Closed pencil point with 0.3mm side port	
M9279	19G/23G	1.1 x 0.64mm Lachrymal Cannula. Straight reducing	
M9278	25G x 7/8"	0.5 x 22mm Lachrymal Cannula. Straight with 0.2mm side port	

DCR Sets

Sterimedix DCR sets offer the latest manufacturing technology and design. Improved attachment of the silicone rod to the malleable probes offers the practitioner the best overall device for DCR surgery.







Code	Gauge and size	Description
M9290	20G x 1 3/4"	0.9 x 45mm Straight DCR Set with 390mm Silicone Rod
M9291	20G x 1 3/4"	0.9 x 45mm Angled DCR Set with 390mm Silicone Rod











VITREORETINAL PRODUCTS

Each instrument is designed for a specific procedure in vitreoretinal surgery.


Backflush and Flute Handles

Code	Gauge, size and Description	Tip
M4489	100mm long 'Rigid' Backflush Handle for Passive Aspiration with finger control and a luer lok cannula connection	
M4491	 'Rigid' Backflush Handle and 20G Silicone Brush set	
M4492	 'Rigid' Backflush Handle and 23G Silicone Brush set	
M4496	 'Rigid' Backflush Handle and 20G Aspiration Cannula set	
M4497	 'Rigid' Backflush Handle and 23G Aspiration Cannula set	














Code	Gauge and size	Description	Tip
M4481	 23G x 1¼" 0.64 x 32mm	Flute Handle with Cannula	
M4480	 20G x 1¼" 0.9 x 32mm	Flute Handle with Cannula	
M4484	 20G x 1¼" 0.9 x 32mm	Active Flute Handle with Cannula	

Code	Gauge, size and Description	Tip
M4493	Charles Flute Handle	
M4494	 Charles Flute Handle with 20G Aspiration Cannula set	
M4495	 Charles Flute Handle with 23G Aspiration Cannula set	





Code	Description
M4498	100mm long heavy liquid infusion handle with finger control and a luer lok cannula connection



Vitreoretinal Cannulae

Code	Gauge and size	Description	Tip
NEW M4470C	 20G/23G 1¼" 0.9/0.64 x 32mm	Curved Reducing Cannula. Reduced 3mm from end of tip	
M4470	 20G/23G 1¼" 0.9/0.64 x 32mm	Reducing Cannula. Reduced 3mm from end of tip	
M4477	 25G x 1¼" 0.5 x 32mm	Aspiration Cannula	
M4478	 23G x 1¼" 0.64 x 32mm	Aspiration Cannula	
M4471	 20G x 1¼" 0.9 x 32mm	Aspiration Cannula	
M4474	 20G/25G 0.9/0.5mm	Perfluorocarbon Coaxial Cannula	
M4499	 23G/30G 0.64/0.3mm	Perfluorocarbon Coaxial Cannula	

Subretinal Fluid Cannulae

Code	Gauge and size	Description	Tip
M4410	 20G 0.9mm	Subretinal Fluid Cannula (Straight). Silicone sleeve extends 6mm from end of cannula	
M4411	 20G 0.9mm	Subretinal Fluid Cannula (Straight). Silicone sleeve extends 3mm from end of cannula	

High Viscosity Injection/Removal

Code	Gauge and size	Description	Tip
M4450	 19G x ½" 1.1 x 12.5mm	Silicone Oil Removal Cannula. Bevelled tip	
M4451	 16G x ¾" 1.6 x 10mm	High Viscosity Injector. Bevelled tip	
M4473	 18G x 1¼" 1.24 x 32mm	Heavy Oil Removal Cannula. Bevelled tip	

High Pressure Infusion Cannulae

Code	Gauge and size	Description
150mm overall length		
M4469	23G 4.0mm	Infusion Cannula
M4460	20G 2.5mm	Infusion Cannula
M4461	20G 4.0mm	Infusion Cannula
M4462	20G 6.0mm	Infusion Cannula


M4463	125mm (5")	High Pressure Extension Tube
-------	------------	------------------------------

Code	Gauge and size	Description
300mm overall length		
M4460A	20G 2.5mm	Infusion Cannula
M4461A	20G 4.0mm	Infusion Cannula
M4462A	20G 6.0mm	Infusion Cannula

Self-retaining Infusion Cannulae

Code	Gauge and size	Description
Silicone Tubing. 300mm long		
M4464	20G 4.3mm	Infusion Cannula (Angled)

M4465	20G 4.3mm	Infusion Cannula (Straight)
M4466	20G 2.5mm	Infusion Cannula (Straight)
M4467	20G 6.0mm	Infusion Cannula (Straight)


Code	Gauge and size	Description	Tip
M4465A	20G 4.3mm	Infusion Cannula (Straight)	
M4466A	20G 2.5mm	Infusion Cannula (Straight)	
M4467A	20G 6.0mm	Infusion Cannula (Straight)	
M4468	23G 4.3mm	Infusion Cannula (Straight)	

Silicone Brush

Code	Gauge and size	Description	Tip	
M4459	● 23G x 1¼"	0.64 x 32mm	Silicone brush. Silicone extends 0.6mm beyond end of cannula	
M4475	● 20G x 1¼"	0.9 x 32mm	Silicone brush. Silicone extends 1.25mm beyond end of cannula	
M4475A	● 20G x 1¼"	0.9 x 32mm	Silicone brush. Silicone extends 0.50mm beyond end of cannula	

Soft Tip Cannulae

Code	Gauge and size	Description	Tip	
M4458	● 23G x 1¼"	0.64 x 32mm	Soft Tip. Silicone extends 2mm beyond end of cannula	




DRAPES & FLUID COLLECTION

Sterimedix drapes are offered in various sizes and materials and are suitable for the left or right eye. They are either waterproof or moisture resistant and the incise offers sufficient adhesive without being too aggressive on removal.

Code	Size & Description
------	--------------------

D1061 122cm x 127cm
(48" x 50")


Incise Drape
with Fluid
Collection Pouch.
Flexible Polymer
Incise Area:
110 x 120mm



Code	Size & Description
------	--------------------


D1063 122cm x 127cm
(48" x 50")

Aperture Drape
with Fluid
Collection Pouch.
Non-woven
Incise Area:
110 x 120mm
Aperture :
44 x 75mm



D1062 122cm x 127cm
(48" x 50")

Incise Drape
with Fluid
Collection Pouch.
Non-woven
Incise Area:
110 x 120mm



PRODUCT CODE	SIZE	DESCRIPTION	PAGE	No. PER BOX
D1061	122 x 127cm	Incise Drape with Pouch	24	5
D1062	122 x 127cm	Incise Drape with Pouch	24	5
D1063	122 x 127cm	Aperture Drape with Pouch	24	5
M0275	25G x 1½"	Standard Retrobulbar Needle	7	10
M0278	19G x 1"	Sub-Tenon Anaesthesia Cannula	7	10
M0279 NEW	20G x 1"	Sub-Tenon Anaesthesia Cannula	7	10
M0285	30G x ¾"	Facial Infiltration Needle	7	10
M0637	25G x 1½"	Atkinson Point Retrobulbar Needle	7	10
M0638	23G x 1½"	Atkinson Point Retrobulbar Needle	7	10
M0641	27G x ¾"	Peribulbar Needle	7	10
M0641A	27G x 1¼"	Peribulbar Needle	7	10
M0642	25G x 1¼"	Peribulbar Needle	7	10
M0642A	25G x ¾"	Peribulbar Needle	7	10
M0643	23G x 1¼"	Peribulbar Needle	7	10
M2270	30G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2270A NEW	30G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2272	26G x ¾"	Visco-Elastic Cannula (Angled)	9	10
M2273	30G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2273A	27G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2273B	20G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2273C	27G x ¾"	Sauter Curved Hydrodissector	12	10
M2273D	25G x ¾"	Sauter Curved Hydrodissector	12	10
M2273E	27G x ¾"	Helsinki Curved Hydrodissector	12	10
M2273F NEW	25G x ¾"	Helsinki Hydrodissector	12	10
M2273G	23G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2273H	27G x ¾"	Curved Hydrodissection Cannula	12	10
M2274	25G x ¾"	Kellan Curved Hydrodissection Cannula	12	10
M2275	25G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2279 NEW	27G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2639	25G x 1"	Rycroft Anterior Chamber Cannula	8	10
M2639A	25G x ¾"	Visco-Elastic Cannula (Angled)	9	10
M2642	21G x ¾"	Rycroft Anterior Chamber Cannula	8	10
M2656	20G x ¾"	Micro Etched A C Maintainer (Dahan)	8	5
M2657	20G x ¾"	Self Retaining A C Maintainer (Lewicky)	8	5
M2658	23G x ¾"	Self Retaining A C Maintainer (Lewicky)	8	5
M2801	25G x ¾"	Lin Lasik Cannula 2 Ports	6	10
M2802	25G x 1½"	Buratto Lasik Irrigating Cannula 2 Ports	6	10
M2804	25G x ¾"	Lasik Irrigating Cannula 6 Ports	6	10
M2817	27G x ¾"	Lasik Irrigating Cannula 4 Ports	6	10
M2819	27G x ¾"	Lasik Irrigating Cannula	6	10
M2822	27G x ¾"	Lasik Irrigating Cannula 4 Ports	6	10
M2823 NEW	25G x ¾"	Lasik Irrigating Cannula (Formed)	6	10
M2824 NEW	23G x ¾"	Lasik Irrigating Cannula (Formed)	6	10
M3000D	25G x ¾"	Capsularhexis Irrigating Cystotome (Formed) (Short curve)	11	10
M3000E	27G x ¾"	Capsularhexis Irrigating Cystotome (Formed) (Short curve)	11	10
M3000F	27G x ¾"	Capsularhexis Irrigating Cystotome (Formed)	11	10
M3000G	25G x ¾"	Capsularhexis Irrigating Cystotome (Formed)	11	10

PRODUCT CODE	SIZE	DESCRIPTION	PAGE	No. PER BOX
M3001	27G x ¾"	Endocapsular Irrigating Cystotome (Joyce) (Curved)	11	10
M3602	25G x ¾"	Side Cutting Cystotome (Formed)	10	10
M3610	25G x ¾"	Irrigating Cystotome (Formed)	10	10
M3610B	27G x ¾"	Irrigating Cystotome (Formed)	10	10
M3610C	25G x ¾"	Irrigating Cystotome (Formed) (45° Angled Tip)	10	10
M3610S	25G x ¾"	Irrigating Cystotome (Formed with Short Tip)	10	10
M3615	27G x ¾"	Irrigating Cystotome (Formed)	10	10
M3615A	25G x ¾"	Irrigating Cystotome (Formed)	10	10
M3626	27G x ¾"	Irrigating Cystotome (Angled)	11	10
M3626A	25G x ¾"	Irrigating Cystotome (Angled)	11	10
M3627	25G x ¾"	Irrigating Cystotome (Formed) (Short Curve)	10	10
M3627A	27G x ¾"	Irrigating Cystotome (Formed) (Short Curve)	10	10
M3629	25G x ¾"	Irrigating Cystotome (Formed, Reverse Tip)	11	10
M3629A	27G x ¾"	Irrigating Cystotome (Formed, Reverse Tip)	11	10
M3630	30G x ¾"	Irrigating Cystotome (Formed)	10	10
M3899	25G x ¾"	Visco-Expression Cannula ('U' Curved with Flat Tip)	9	10
M3899A	25G x ¾"	Visco-Expression Cannula ('U' Curved with Round Tip)	9	10
M3900	25G x ¾"	Hydrodissection Cannula (Angled)	12	10
M3901	25G x ¾"	Hydrodissection Cannula ('J' Shaped, Flattened Tip)	13	10
M3902	27G x ¾"	Hydrodissection Cannula (Angled)	12	10
M3904	27G x ¾"	Chang Hydrodissection Cannula	12	10
M3906A	25G x ¾"	Micro J Shaped Hydrodissector	13	10
M3907	25G x 7/8"	Tapered Hydrodelineator	9	10
M3908	25G x 7/8"	Cortical Cleaving Cannula	9	10
M3909	23G x ¾"	Visco-Elastic Cannula (Angled)	9	10
M3910	21G x ¾"	Visco-Elastic Cannula (Straight)	9	10
M3910A	21G x ¾"	Visco-Elastic Cannula (Angled)	9	10
M3911	27G x ¾"	Cortical Cleaving Hydrodissector (Angled, Flattened Tip)	13	10
M3911A	27G x ¾"	Cortical Cleaving Hydrodissector (Angled)	13	10
M3912 NEW	23G x ¾"	Mackool Hydrodissector	12	10
M3913 NEW	25G x ¾"	Daya Hydrodissector	12	10
M3914 NEW	23G x ¾"	Daya Hydrodissector	12	10
M3917 NEW	27G x ¾"	Visco-Elastic Cannula (Formed)	9	10
M3926A	27G x ¾"	Visco-Elastic Cannula (Angled)	9	10
M3927 NEW	27G x ¾"	Visco-Elastic Cannula (Angled)	9	10
M3992	30G x ¾"	Nucleus Hydrodissector	12	10
M4410	20G	Subretinal Fluid Cannula 6mm (Straight)	22	10
M4411	20G	Subretinal Fluid Cannula 3mm (Straight)	22	10
M4450	19G x ½"	Silicone Oil Removal Cannula	22	10
M4451	16G x ¾"	High Viscosity Injector	22	10
M4458	23G x 1¼"	Soft Tip Cannula	23	10
M4459	23G x 1¼"	Silicone Brush	23	10
M4460	20G x 2.5mm	Infusion Cannula x 150mm	23	10
M4460A	20G x 2.5mm	Infusion Cannula x 300mm	23	10

PRODUCT CODE	SIZE	DESCRIPTION	PAGE	No. PER BOX
M4461	20G x 4.0mm	Infusion Cannula x 150mm	23	10
M4461A	20G x 4.0mm	Infusion Cannula x 300mm	23	10
M4462	20G x 6.0mm	Infusion Cannula x 150mm	23	10
M4462A	20G x 6.0mm	Infusion Cannula x 300mm	23	10
M4463	20G x 125mm	High Pressure Extension Tube	23	10
M4464	20G x 4.3mm	Self Retaining Angled Infusion Cannula x 300mm	23	10
M4465	20G x 4.3mm	Self Retaining Infusion Cannula x 300mm	23	10
M4465A	20G x 4.3mm	Self Retaining Infusion Cannula x 150mm	23	10
M4466	20G x 2.5mm	Self Retaining Infusion Cannula x 300mm	23	10
M4466A	20G x 2.5mm	Self Retaining Infusion Cannula x 150mm	23	10
M4467	20G x 6.0mm	Self Retaining Infusion Cannula x 300mm	23	10
M4467A	20G x 6.0mm	Self Retaining Infusion Cannula x 150mm	23	10
M4468	23G x 4.3mm	Self Retaining Infusion Cannula (Straight) x 150mm	23	10
M4469	23G x 4.0mm	Infusion Cannula (Straight) x 150mm	23	10
M4470	20G/23G	Reducing Cannula	22	10
M4470C NEW	20G/23G	Curved Reducing Cannula	22	10
M4471	20G x 1½"	Aspiration Cannula	22	10
M4473	18G x 1½"	Heavy Oil Removal Cannula	22	10
M4474	20G/25G	Perfluorocarbon Cannula	22	10
M4475	20G x 1½"	Silicone Brush	23	10
M4475A	20G x 1½"	Silicone Brush	23	10
M4477	25G x 1½"	Aspiration Cannula	22	10
M4478	23G x 1½"	Aspiration Cannula	22	10
M4480	20G x 1½"	Flute Handle with Cannula	21	5
M4481	23G x 1½"	Flute Handle with Cannula	21	5
M4484	20G x 1½"	Active Flute Handle with Cannula	21	5
M4489	100mm	'Rigid' Backflush Handle	21	5
M4491	100mm 20G x 1½"	'Rigid' Backflush Handle + 20G Silicone Brush	21	5
M4492	100mm 23G x 1½"	'Rigid' Backflush Handle + 23G Silicone Brush	21	5
M4493	110mm	Charles Flute Handle	21	5
M4494	110mm 20G x 1½"	Charles Flute Handle + 20G Aspiration Cannula	21	5
M4495	110mm 23G x 1½"	Charles Flute Handle + 23G Aspiration Cannula	21	5
M4496	100mm 20G x 1½"	'Rigid' Backflush Handle + 20G Aspiration Cannula	21	5
M4497	100mm 23G x 1½"	'Rigid' Backflush Handle + 23G Aspiration Cannula	21	5
M4498	110mm	Heavy Liquid Infusion Cannula	21	5
M4499	23G/30G	Perfluorocarbon Cannula	22	10
M4619	25G x 1½"	Irrigating Vectis (Single Port)	14	10
M4620	25G x 1½"	Irrigating Vectis (Double Port)	14	10
M4621	27G x 1½"	Small Incision Vectis (Single Port)	14	10






PRODUCT CODE	SIZE	DESCRIPTION	PAGE	No. PER BOX
M5405	21G x ¾"	Bimanual Irrigation Cannula	17	10
M5409	21G x ¾"	Bimanual Irrigation Cannula (2 x 0.5mm Ports)	17	10
M5414	21G x ¾"	Bimanual Irrigation Cannula (2 x 0.5mm Ports)	17	10
M5501	23G x ¾"	Aspiration Handle 0.4mm Side Port and Micro Etched Tip (Curved)	16	5
M5502	23G x ¾"	Aspiration Handle 0.3mm Side Port and Micro Etched Tip (Curved)	16	5
M5505	21G x ¾"	Irrigation Handle (Curved)	16	5
M5509	21G x ¾"	Irrigation Handle with 2 x 0.5mm Ports (Straight)	16	5
M5510	20G x ¾"	Irrigation Handle with 2 x 0.5mm Ports (Curved)	16	5
M5511	20G x ¾"	Aspiration Handle 0.3mm Side Port and Micro Etched Tip (Curved)	16	5
M5512	21G x ¾"	Aspiration Handle 0.3mm Side Port and Micro Etched Tip (Curved)	16	5
M5513 NEW	23G x ¾"	Irrigation Handle with 2 x 0.4mm Ports (Curved)	16	5
M5514	21G x ¾"	Irrigation Handle with 2 x 0.5mm Ports (Curved)	16	5
M5515	21G x ¾"	Aspiration Handle with 0.3mm side Port (Curved)	16	5
M5518	23G x ¾"	Aspiration Handle 0.3mm Side Port (Curved)	16	5
M5519	23G x ¾"	Irrigation Handle Ultra Thin Wall with Bevelled Tip (Curved)	16	5
M5520	21G x ¾"	Aspiration Handle Thin Wall with 0.4mm Side Port and Micro Etched Tip (Curved)	16	5
M5521	21G x ¾"	Irrigating Handle Thin Wall with 2 x 0.6mm Side Ports	16	5
M5522	21G x ¾"	Irrigating Handle Thin Wall with Open End	16	5
M5523	23G x ¾"	Irrigating Handle Ultra Thin Wall with Open End (Curved)	16	5
M5531	21G x ¾"	Bimanual I/A System Combining: Irrigating Handle Thin Wall with 2 x 0.6mm Side Ports	17	5 Pairs
	21G x ⅝"	Aspiration Handle Thin Wall with 0.4mm Side Port and Micro Etched Tip (Curved)	17	
M5532	21G x ¾"	Bimanual I/A System Combining: Irrigating Handle Thin Wall with Open End	17	5 Pairs
	21G x ¾"	Aspiration Handle Thin Wall with 0.4mm Side Port and Micro Etched Tip (Curved)	17	
M5535 NEW	23G x ¾"	Bimanual I/A System Combining: Irrigating Handle Thin Wall with 2 x 0.4mm Side Ports	17	5 Pairs
	23G x ¾"	Aspiration Handle Thin Wall with 0.3mm Side Port and Micro Etched Tip (Curved)	17	
M5580	16G/21G	I/A Handpiece (Straight Tip)	18	5
M5581	16G/21G	I/A Handpiece (Curved Tip)	18	5
M5582	16G/21G	I/A Handpiece (45° Angled Tip)	18	5
M5583	16G/21G	I/A Handpiece (60° Angled Tip)	18	5
M5584	16G/21G	I/A Handpiece (90° Angled Tip)	18	5
M5620	23G x 1"	Binkhorst Hook (Right)	14	10
M5621	23G x 1"	Binkhorst Hook (Left)	14	10
M5622	25G x ¾"	Binkhorst Hook (Right)	14	10
M5623	25G x ¾"	Binkhorst Hook (Left)	14	10

PRODUCT CODE	SIZE	DESCRIPTION	PAGE	No. PER BOX
M6524	21G x 5/8"	Bimanual Aspirating Cannula	17	10
M6525	21G x 5/8"	Cortex Extractor 0.3mm Port	14	10
M6575	23G x 5/8"	Cortex Extractor 0.3mm Port	14	10
M6575A	23G x 5/8"	Cortex Extractor 0.4mm Port	14	10
M6576	23G x 5/8"	Cortex Extractor 0.3mm Port	14	10
M6576A	23G x 5/8"	Cortex Extractor 0.3mm Port	14	10
M6594	23G x 5/8"	Cortex Extractor (Charleux)	14	10
M6595	23G x 5/8"	Cortex Extractor (Charleux)	14	10
M6596	24G x 1/2"	Cortex Extractor (Charleux)	14	10
M6600	23G/23G	Standard Simcoe (Aspiration through 0.4mm Side Port)	15	5
M6601	23G/23G	Reverse Simcoe (Aspiration through 0.4mm Side Port)	15	5
M6605	23G/23G	Standard Simcoe (Aspiration through 0.3mm Side Port)	15	5
M6611	21G x 1 1/2"	Single Cannula (Straight)	15	10
M6634	18G/21G	Coaxial I/A Cannula (Straight)	15	5
M6634B	18G/21G	Coaxial I/A Cannula (Angled)	15	5
M7500	3.5mm/4.0mm	Scleral Marker	6	10
M8603	21G x 5/8"	Capsule Polisher	19	10
M8604	21G x 1"	Posterior Capsule Polisher	19	10
M8604A	23G x 5/8"	Posterior Capsule Polisher	19	10
M8604B	21G x 5/8"	Posterior Capsule Polisher	19	10
M8605	27G x 5/8"	Silicone Tip Polisher / Manipulator	19	10
M8605A	23G x 5/8"	Silicone Tip Polisher / Manipulator	19	10
M8608	27G x 5/8"	Flexible Tipped Polisher	19	10
M8620	23G x 1"	Olive Tipped Cannula	19	10
M8620C	23G x 1"	Olive Tipped Cannula (Curved)	19	10
M8621	23G x 1"	Olive Tipped Cannula (Curved) with Micro-etched Olive	19	10
M8622	25G x 5/8"	Olive Tipped Cannula (Angled) with Micro-etched Olive	19	10
M8645	23G x 5/8"	Posterior Scraper	19	10
M9272	23G x 1/2"	Lachrymal Cannula (Short Straight)	20	10
M9274	21G x 1 1/2"	Lachrymal Cannula (Curved)	20	10
M9275	26G x 3/4"	Lachrymal Cannula (Straight)	20	10
M9276	26G x 1 1/4"	Lachrymal Cannula (Curved)	20	10
M9277	25G x 1 1/8"	Lachrymal Cannula (Curved)	20	10
M9278	25G x 5/8"	Lachrymal Cannula (Straight)	20	10
M9279	19G/23G	Lachrymal Cannula / Dilator	20	10
M9280	19G x 1"	Washout Cannula (Bishop Harmon)	8	10
M9281	16G x 1"	Washout Cannula	8	10
M9290	20G x 1 1/4"	Straight DCR Set with 390mm Silicone Rod	20	5
M9291	16G x 1"	Angled DCR Set with 390mm Silicone Rod	20	5

Product key

Colour

Gauge

Yellow			30G
Grey			27G
Brown			26G
Orange			25G
Purple			24G
Blue			23G
Black			22G
Green			21G
Yellow			20G
Cream			19G
Pink			18G
White			16G

All hub colours comply with ISO6009

KIT PACKERS

All Sterimedix products are available to kit packers non-sterile.

Drapes are provided in bulk, non-sterile

For further information, please contact us at:


sales@sterimedix.com





STERIMEDIX  [®]

delivering **innovation** in ophthalmic surgery



Sterimedix Limited 1 Madeley Road North Moons Moat Redditch Worcestershire B98 9NB UK
T : +44 (0)1527 501480 F : +44 (0)1527 501491 e : sales@sterimedix.com

www.sterimedix.com

Sterimedix policy is one of continuous improvement and we therefore reserve the right to alter specifications and prices without prior notice.
All products are sold subject to the company's conditions of sale.

STERIMEDIX  [®] is a registered trademark of Sterimedix Limited

08.2011