









ANAESTHESIA



Retrobulbar

For injecting anaesthetic into the muscle cone.

Code	Gauge and size	Description	Tip
M0275	 25G x 1½"	0.5 x 38mm Retrobulbar Tri-Facet	
M0285	 30G x 5⁄8"	0.3 x 16mm Facial Infiltration	
M0637	 25G x 1½"	0.5 x 38mm Retrobulbar (Atkinson)	
M0638	 23G x 1½"	0.64 x 38mm Retrobulbar Thin Wall (Atkinson)	






Sub Tenon Cannula

For injecting anaesthetic into the sub tenon space.

Code	Gauge and size	Description	Tip
M0278	 19G x 1"	1.1 x 25mm Sub Tenon Cannula Curved with Flattened Tip	

Peribulbar

For injecting anaesthetic around the globe.

Code	Gauge and size	Description	Tip
M0641A	 27G x 1¼"	0.4 x 32mm Peribulbar	
M0642	 25G x 1¼"	0.5 x 32mm Peribulbar	
M0643	 23G x 1¼"	0.64 x 32mm Peribulbar	
M0641	 27G x 7⁄8"	0.4 x 22mm Peribulbar	
M0642A	 25G x 7⁄8"	0.5 x 22mm Peribulbar	