














Refractive Cannulas

Refractive Cannulas

Products designed for laser refractive procedures.

Code	Gauge and size	Description	Tip
M2804	 25G x 7/8" 0.5 x 22mm	Irrigating Lasik Cannula. Formed 7mm from tip. Six 0.15mm ports at 2, 4 and 6mm from tip	
M2817	 27G x 7/8" 0.4 x 22mm	Irrigating Lasik Cannula. Formed 9mm from tip. Four 0.15mm ports at 2 and 4mm from tip	
M2822	 25G x 7/8" 0.5 x 22mm	Irrigating Lasik Cannula. Formed 7mm from tip. Six 0.15mm ports at 2, 4 and 6mm from closed tip	
 M2824	 26g x 7/8" 0.45 x 22mm	Irrigating Lasik Cannula. Formed 9.5mm from spatulated tip (Slade)	
M2801	 25G x 7/8" 0.5 x 22mm	Irrigating Lasik Cannula. Flattened tip. Angled with two 0.25mm ports (Lin)	
M2802	 25G x 1 1/8" 0.5 x 29mm	Irrigating Lasik Cannula. Formed with two 0.20mm side ports (Buratto)	

Corneal Surgery Cannulas

Code	Gauge and size	Description	Tip
 M8201	 27G x 7/8" 0.4 x 22mm	DALK Cannula. Formed 3.5mm from tip with 0.2mm port on under-side	