











Refractive Cannulas / Scleral Markers

Refractive Cannulas

Products designed for irrigating the cornea and stroma during laser refractive procedures.

| Code | Gauge and size | Description | Tip |
|-------|---|---|---|
| M2804 |  25G x 7/8" 0.5 x 22mm | Irrigating Lasik Cannula. Formed 7mm from tip. Six 0.15mm ports at 2, 4 and 6mm from tip |  |
| M2817 |  27G x 7/8" 0.4 x 22mm | Irrigating Lasik Cannula. Formed 9mm from tip. Four 0.15mm ports at 2 and 4mm from tip |  |
| M2822 |  25G x 7/8" 0.5 x 22mm | Irrigating Lasik Cannula. Formed 7mm from tip. Six 0.15mm ports at 2, 4 and 6mm from closed tip |  |
| M2824 |  26g x 7/8" 0.45 x 22mm | Irrigating Lasik Cannula. Formed 9.5mm from spatulated tip (Slade) |  |
| M2801 |  25G x 7/8" 0.5 x 22mm | Irrigating Lasik Cannula. Flattened tip. Angled with two 0.25mm ports (Lin) |  |
| M2802 |  25G x 1 1/8" 0.5 x 29mm | Irrigating Lasik Cannula. Formed with two 0.2mm side ports (Buratto) |  |
| M3914 |  23G x 7/8" 0.64 x 22mm | Irrigating Lasik Cannula. Formed 10.5mm from flattened tip. |  |