



Anaesthetic Needles and Cannulas



Anaesthetic Cannulas

For injecting anaesthetic into the muscle cone.

Code	Gauge and size	Description	Tip
			
M0278	● 19G x 1" 1.1 x 25mm	Sub Tenon Cannula. Curved with flattened tip	
M0279	● 20G x 1" 0.9 x 25mm	Sub Tenon Cannula. Curved with flattened tip	

Anaesthetic Needles

For injecting anaesthetic into the sub tenon space and around the globe.

Code	Gauge and size	Description	Tip
			
M0275	● 25G x 1½" 0.5 x 38mm	Retrobulbar. Tri-facet	
M0637	● 25G x 1½" 0.5 x 38mm	Retrobulbar (Atkinson)	
M0638	● 23G x 1½" 0.64 x 38mm	Retrobulbar (Atkinson)	
			
M0641A	● 27G x 1¼" 0.4 x 32mm	Peribulbar (Atkinson)	
M0642	● 25G x 1¼" 0.5 x 32mm	Peribulbar (Atkinson)	
M0643	● 23G x 1¼" 0.64 x 32mm	Peribulbar (Atkinson)	
M0641	● 27G x 7⁄8" 0.4 x 22mm	Peribulbar (Atkinson)	
M0642A	● 25G x 7⁄8" 0.5 x 22mm	Peribulbar (Atkinson)	