

Anterior Chamber Maintenance Cannulas

Anterior Chamber Maintainers


Self-retaining cannulas designed for the continual infusion of fluid to maintain the AC during cataract surgery or anterior vitrectomy.

Code	Gauge and size	Description	Tip
			
M2658	23G x 9/64"	0.64 x 3.5mm	Anterior Chamber Maintainer. 150mm silicone tubing. Annular Rings allow cannula to be self retaining (Lewicky)
M2657	20G x 9/64"	0.9 x 3.5mm	Anterior Chamber Maintainer. 150mm silicone tubing. Annular Rings allow cannula to be self retaining (Lewicky)



Anterior Chamber Maintenance Cannulas


For maintaining and forming the anterior chamber during surgery by injecting or removing viscoelastic materials, BSS or air.

Code	Gauge and size	Description	Tip
			
M2639	25G x 1"	0.5 x 25mm	Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)
M2270	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 4mm from tip (Rycroft)
M2270A	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 8mm from tip (Rycroft)
M2273	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)
M2273A	27G x 7/8"	0.4 x 22mm	Anterior Chamber Cannula. Angled 35°, 4.5mm from tip (Rycroft)
M2279	27G x 7/8"	0.4 x 22mm	Anterior Chamber Cannula. Angled 45°, 11mm from tip (Rycroft)
M2275	25G x 7/8"	0.5 x 22mm	Anterior Chamber Cannula. Angled 35°, 5mm from tip (Rycroft)
M2273G	23G x 7/8"	0.64 x 22mm	Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)
M2642	21G x 7/8"	0.81 x 22mm	Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)
M2273B	20G x 7/8"	0.9 x 22mm	Anterior Chamber Cannula. Angled 45°, 7.5mm from tip (Rycroft)



Transzonular Injection Cannula

For the intraoperative injection of steroid-antibiotic combinations.

Code	Gauge and size	Description	Tip
			
M2280	27G x 7/8"	0.4 x 22mm	Transzonular Injection Cannula. Angled 33°, 2mm from tip.

