







Viscoelastic Cannulas

For injection of viscoelastic products into the anterior or posterior chamber.

Code	Gauge and size	Description	Tip	
	M3926A	27G x 7/8" 0.4 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip	
M3927	27G x 7/8" 0.4 x 22mm	Viscoelastic Cannula. Angled at 25°, 8mm from tip		
M2272	26G x 7/8" 0.45 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip		
M2639A	25G x 7/8" 0.5 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip		
M3909	23G x 7/8" 0.64 x 22mm	Viscoelastic Cannula. Angled at 45°, 9mm from tip		
M3910A	21G x 7/8" 0.81 x 22mm	Viscoelastic Cannula. Angled at 45°, 9mm from tip		
	 M3910	21G x 7/8" 0.81 x 22mm		Viscoelastic Cannula. Straight
	 M3917	27G x 7/8" 0.4 x 22mm		Viscoelastic Cannula. Angled and formed 11mm from tip