










Capsule Polishing Cannulas

Capsule Polishing Cannulas









Designed for polishing and cleaning the residual cortex from the capsular bag. Silicone tipped cannulas provide gentle polishing and atraumatic insertion and removal. Micro-etched cannulas are slightly roughened for gentle polishing.

Code	Gauge and size	Description	Tip
	 21G x 1" 0.81 x 25mm	Capsule Polisher. Angled 12.5mm from tip. Aspiration through 0.3mm side port. Micro-etched tip	

Silicone Tipped Polishers

Code	Gauge and size	Description	Tip
	 27G x 7/8" 0.4 x 22mm	Silicone tipped polisher. Angled	
	 23G x 7/8" 0.64 x 22mm	Silicone tipped polisher. Angled	
	 27G x 7/8" 0.4 x 22mm	Silicone tipped polisher. Curved through 45°	

Olive Tipped Cannulas

Code	Gauge and size	Description	Tip
	 23G x 1" 0.64 x 25mm	Olive tipped polisher. Straight	
	 23G x 1" 0.64 x 25mm	Olive tipped polisher. Curved	
	 23G x 1" 0.64 x 25mm	Olive tipped polisher. Curved. Micro-etched tip	