











































Irrigating/Aspirating Cannulas and Handpieces

Coaxial Cannulas and Handpieces

Code	Gauge and size	Description	Tip	
	M6634	18G/21G 1.24/0.81mm	Coaxial I & A Cannula. 0.3mm side aspiration port on inner cannula	
	M6611	 21G x 1½" 0.81 x 38mm	I/A Straight Cannula	
	 M5590	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. Straight tip	
	 M5592	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 45° angled tip	
	 M5593	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 60° angled tip	
	 M5594	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 90° angled tip	
	 M5591	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. Curved tip	
	M5580	16/21G 1.6/0.81mm	I/A Handpiece. Straight tip	
	M5582	16/21G 1.6/0.81mm	I/A Handpiece. 45° angled tip	
	M5583	16/21G 1.6/0.81mm	I/A Handpiece. 60° angled tip	
	M5584	16/21G 1.6/0.81mm	I/A Handpiece. 90° angled tip	
	M5581	16/21G 1.6/0.81mm	I/A Handpiece. Curved tip.	
	 M5560	21G	Soft Outer I/A Handpiece. Straight tip	
	 M5561	21G	Soft Outer I/A Handpiece. 45° angled tip	
	 M5562	21G	Soft Outer I/A Handpiece. 60° angled tip	
	 M5563	21G	Soft Outer I/A Handpiece. 90° angled tip	

Irrigating/Aspirating Cannulas and Handpieces



new

M5570

24G

Soft Outer I/A Handpiece for 2.2mm incisions. Straight tip



new

M5571

24G

Soft Outer I/A Handpiece for 2.2mm incisions. 45° angled tip



new

M5572

24G

Soft Outer I/A Handpiece for 2.2mm incisions. 60° angled tip



new

M5573







24G

Soft Outer I/A Handpiece for 2.2mm incisions. 90° angled tip



Simcoe Irrigating / Aspirating Cannulas

For removing cortical debris or viscoelastic substances while maintaining the anterior chamber during extracapsular and intercapsular/endocapsular cataract procedures.

Code	Gauge and size	Description	Tip	
	M6600	23G/23G 0.64/0.64mm	I/A Cannula Standard. 0.4mm port. 200mm silicone tubing 2.65 OD x 0.65 ID (Simcoe)	
	M6601	23G/23G 0.64/0.64mm	I/A Cannula Reverse. 0.4mm port. 200mm silicone tubing 2.65 OD x 0.65 ID (Simcoe)	
	M6605	23G/23G 0.64/0.64mm	I/A Cannula Standard. 0.3mm port. 200mm silicone tubing 2.00 OD x 1.00 ID (Simcoe)	

Irrigating Microsurgical Handle

Code	Gauge and size	Description	Tip
	M1000	Irrigating Microsurgical Handle with male luer lock connector	

new