










Capsule Polishing Cannulas

Capsule Polishing Cannulas






Designed for polishing and cleaning the residual cortex from the capsular bag. Silicone tipped cannulas provide gentle polishing and atraumatic insertion and removal. Micro-etched cannulas are slightly roughened for gentle polishing. Olive tipped cannulas allow polishing over a greater surface area with a reduced risk of posterior capsule tear.

Code	Gauge and size	Description	Tip
			
M8604A	23G x 7/8" 0.64 x 22mm	Capsule Polisher. Angled 9mm from tip. Aspiration through 0.3mm side port. Micro-etched tip	
M8604	21G x 1" 0.81 x 25mm	Capsule Polisher. Angled 12.5mm from tip. Aspiration through 0.3mm side port. Micro-etched tip	
			
 M8604B	21G x 7/8" 0.81 x 22mm	Capsule Polisher. Angled 9mm from tip. Aspiration through 0.3mm side port. Micro-etched tip.	
 M8645	23G x 7/8" 0.64 x 22mm	Capsule Scraper. Angled 0.5 & 7mm from tip. Micro-etched distal end.	

Silicone Tipped Polishers

Code	Gauge and size	Description	Tip
			
M8605	27G x 7/8" 0.4 x 22mm	Silicone tipped polisher. Angled	
M8605A	23G x 7/8" 0.64 x 22mm	Silicone tipped polisher. Angled	
			
M8608	27G x 7/8" 0.4 x 22mm	Silicone tipped polisher. Curved through 45°	

Olive Tipped Cannulas

Code	Gauge and size	Description	Tip
			
 M8622	25G x 1" 0.5 x 25mm	Olive tipped polisher. Angled 10mm from tip. Micro-etched tip.	
			
M8620C	23G x 1" 0.64 x 25mm	Olive tipped polisher. Curved	
M8621	23G x 1" 0.64 x 25mm	Olive tipped polisher. Curved. Micro-etched tip	