

























Irrigating/Aspirating Cannulas and Handpieces

Coaxial Cannulas and Handpieces

For the aspiration of cortex while simultaneously maintaining the AC with a continuous irrigating flow.

Code	Gauge and size	Description	Tip
 M6634	18G/21G 1.24/0.81mm	Coaxial I & A Cannula. 0.3mm side aspiration port on inner cannula	
 M6611	 21G x 1½" 0.81 x 38mm	I/A Straight Cannula	
 M5590	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. Straight tip	
M5592	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 45° angled tip	
M5593	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 60° angled tip	
M5594	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 90° angled tip	
 M5591	17/24G 1.47/0.56mm	I/A Handpiece for 2.2mm incisions. Curved tip	
 M5580	16/21G 1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. Straight tip	
M5582	16/21G 1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. 45° angled tip	
M5583	16/21G 1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. 60° angled tip	
M5584	16/21G 1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. 90° angled tip	
 M5581	16/21G 1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. Curved tip.	
 M5560	21G	Soft Outer I/A Handpiece for 2.7/2.8mm incisions. Straight tip	
M5561	21G	Soft Outer I/A Handpiece for 2.7/2.8mm incisions. 45° angled tip	
M5562	21G	Soft Outer I/A Handpiece for 2.7/2.8mm incisions. 60° angled tip	
M5563	21G	Soft Outer I/A Handpiece for 2.7/2.8mm incisions. 90° angled tip	

Irrigating/Aspirating Cannulas and Handpieces



M5570	24G	Soft Outer I/A Handpiece for 2.2mm incisions. Straight tip	
M5571	24G	Soft Outer I/A Handpiece for 2.2mm incisions. 45° angled tip	
M5572	24G	Soft Outer I/A Handpiece for 2.2mm incisions. 60° angled tip	
M5573	24G	Soft Outer I/A Handpiece for 2.2mm incisions. 90° angled tip	

Simcoe Irrigating / Aspirating Cannulas

For removing cortical debris or viscoelastic substances while maintaining the anterior chamber during extracapsular and intercapsular/endocapsular cataract procedures.

Code	Gauge and size	Description	Tip
------	----------------	-------------	-----



M6600	● 23G/23G	0.64/0.64mm	I/A Cannula Standard. 0.4mm port. 200mm silicone tubing 2.65 OD x 0.65 ID (Simcoe)	
M6602	● 23G/23G	0.64/0.64mm	I/A Cannula Standard. 0.4mm port. 270mm silicone tubing 2.65 OD x 0.65 ID (Simcoe)	



M6605	● 23G/23G	0.64/0.64mm	I/A Cannula Standard. 0.3mm port. 200mm silicone tubing 2 OD x 1 ID (Simcoe)	
-------	-----------	-------------	--	--

Irrigating Microsurgical Handle

Code	Gauge and size	Description	Tip
------	----------------	-------------	-----



M1000		Irrigating Microsurgical Handle with female luer taper and a male luer lock connector	
-------	--	---	--