















# Lachrymal Cannulas and Dacryocystorhinostomy (DCR) Sets


## Lachrymal Cannulas

Cannulas for probing and irrigating the lachrymal ducts

Code	Gauge and size	Description	Tip
 M9275	 26G x 3/4"	0.45 x 20mm Lachrymal Cannula. Short straight	
 M9276	 26G x 1 1/4"	0.45 x 32mm Lachrymal Cannula. Curved	
M9277	 25G x 1 1/8"	0.5 x 29mm Lachrymal Cannula. Curved	
M9274	 21G x 1 1/2"	0.81 x 38mm Lachrymal Cannula. Curved	
 M9279	 19G/23G	1.1 x 0.64mm Lachrymal Cannula. Straight reducing	
 M9278	 25G x 7/8"	0.5 x 22mm Lachrymal Cannula. Straight with 0.2mm side port	

## Dacryocystorhinostomy (DCR) Sets

Sterimedix DCR sets offer the latest manufacturing technology and design. Improved attachment of the silicone rod to the malleable probes offers the practitioner the best overall device for DCR surgery.

Code	Gauge and size	Description	Tip
M9290	20G x 1 3/4"	0.9 x 45mm Straight DCR Set with 390mm silicone rod	
M9291	20G x 1 3/4"	0.9 x 45mm Angled DCR Set with 390mm silicone rod	

