






















Irrigating/Aspirating Cannulas and Handpieces







Bimanual I/A Handpieces

All Bimanual handles have a coloured indicator to show the gauge size in accordance with ISO 6009.

Code	Gauge and size	Description	Tip
			
M5501	 23G x 5/8" 0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.4mm side port and micro-etched tip	
M5502	 23G x 5/8" 0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
M5512	 21G x 5/8" 0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
M5520	 21G x 5/8" 0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.4mm side port and micro-etched tip	
M5511	 20G x 3/4" 0.9 x 19mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
			
M5518	 23G x 5/8" 0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	
M5515	 21G x 7/8" 0.81 x 22mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	
M5519	 23G x 5/8" 0.64 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with 45° bevelled tip	
M5505	 21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Handpiece. Curved with 45° bevelled tip	
M5513	 23G x 5/8" 0.64 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.4mm side ports	
M5521	 21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.5mm side ports	
M5514	 21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Handpiece. Curved with twin 0.5mm side ports	
M5510	 20G x 3/4" 0.9 x 19mm	Bimanual Irrigation Handpiece. Curved with twin 0.5mm side ports	
M5523	 23G x 5/8" 0.64 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with open end	

















Bimanual I/A Systems

A combination set of one Irrigating Handle and one Aspirating Handle manufactured using thin wall stainless steel tubing.

Code	Gauge and size	Description	Tip
			
 M5537 - set	 20G x 3/4" 0.9 x 19mm	Bimanual Irrigation Handpiece. Curved with twin 0.5mm side ports	
	 21G x 7/8" 0.81 x 22mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	




Irrigating/Aspirating Cannulas and Handpieces

Bimanual I/A Systems

Code	Gauge and size	Description	Tip
			
M5531 - set	● 21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.5mm side ports	
	● 21G x 5/8" 0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.4mm side port and micro-etched tip	
			
 M5531 - set	● 21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Handpiece. Curved with 2 x 0.4mm irrigation ports and tapered tip	
	● 21G x 5/8" 0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm aspiration port and micro-etched tapered tip	
			
M5535 - set	● 23G x 5/8" 0.64 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.4mm side ports	
	● 23G x 5/8" 0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
			
M5536 - set	● 23G x 5/8" 0.64 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with 45° bevelled tip	
	● 23G x 5/8" 0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	
			
M5538 - set	● 25G x 5/8" 0.5 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with 45° bevelled tip	
	● 25G x 5/8" 0.5 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	

For more information of other systems available, please contact Sterimedix or your local distributor.

Bimanual Cannulas

Code	Gauge and size	Description	Tip
			
M6524	● 21G x 5/8" 0.81 x 16mm	Bimanual Aspiration Cannula. Curved with 0.3mm top port	
			
M5405	● 21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Cannula. Curved with 45° bevelled tip	