























# Vitreoretinal Cannulas

## Vitreoretinal Cannulas







Code	Gauge and size	Description	Tip
M4470C	 20G/23G x 1¼" 0.9/0.64 x 32mm	Reducing Cannula. Curved. Reduced 3mm from tip	
M4470	 20G/23G x 1¼" 0.9/0.64 x 32mm	Reducing Cannula. Reduced 3mm from tip	
M4477	 25G x 1¼" 0.5 x 32mm	Aspiration Cannula	
M4478	 23G x 1¼" 0.64 x 32mm	Aspiration Cannula	
 M4478B	 23G UTW x 1¼" 0.64 x 32mm	Aspiration Cannula. Ultra Thin Wall	
 M4478C	 23G UTW x 1" 0.64 x 25mm	Aspiration Cannula. Ultra Thin Wall	
M4471	 20G x 1¼" 0.9 x 32mm	Aspiration Cannula	
M4410	 20G 0.9mm	Subretinal Fluid Cannula. Straight. Silicone sleeve extends 6mm from end of cannula	



## Heavy Liquid Infusion Aspiration Cannulas






Code	Gauge and size	Description	Tip
 M4454	 20G x ¼" 0.90 x 7mm	Viscous Fluid Injection Cannula	
 M4453	 25G x ¼" 0.5 x 7mm	Viscous Fluid Injection Cannula	
 M4452	 23G x ¼" 0.64 x 7mm	Viscous Fluid Injection Cannula. Bevelled tip	
M4450	 19G x ½" 1.1 x 12.5mm	Heavy Liquid Infusion Aspiration Cannula. Bevelled tip	

## Dual-Bore Cannulas

Code	Gauge and size	Description	Tip
M4474	 20G/25G 0.9/0.5mm	Perfluorocarbon Coaxial (Dual Bore) Cannula	
M4499	 23G/30G 0.64/0.3mm	Perfluorocarbon Coaxial (Dual Bore) Cannula	
 M4499A	 25G/30G 0.5/0.3mm	Perfluorocarbon Coaxial (Dual Bore) Cannula	

# Vitreoretinal Cannulas

## Silicone Brush Cannulas

Code	Gauge and size	Description	Tip
 M4449	 25G x 1¼" 0.5 x 32mm	Silicone Brush. Silicone extends 0.7mm beyond end of cannula	
M4459	 23G x 1¼" 0.64 x 32mm	Silicone Brush. Silicone extends 0.6mm beyond end of cannula	
M4475	 20G x 1¼" 0.9 x 32mm	Silicone Brush. Silicone extends 1.25mm beyond end of cannula	
M4475A	 20G x 1¼" 0.9 x 32mm	Silicone Brush. Silicone extends 0.60mm beyond end of cannula	

## Soft Tip Cannulas

Code	Gauge and size	Description	Tip
 M4456	 25G x 1¼" 0.5 x 32mm	Soft Tip. Silicone extends 1mm beyond end of cannula	
M4458	 23G x 1¼" 0.64 x 32mm	Soft Tip. Silicone extends 2mm beyond end of cannula	
 M4457	 20G x 1¼" 0.64 x 32mm	Soft Tip. Silicone extends 1mm beyond end of cannula	