

# Vitreoretinal Extension Tube & Self Retaining Infusion Cannulas

## High Pressure Extension Tube

Code

Gauge and size

Description



M4463

125mm

High Pressure Extension Tube



## Self Retaining Infusion Cannulas

For simple insertion and removal without the need for a retaining suture. Designed for use directly into the eye without the need for a trocar.

Code

Gauge and size

Description

Tip



300mm overall length

M4464

20G

4.3mm

Self Retaining Infusion Cannula. Angled



150mm overall length

M4467A

20G

6mm

Self Retaining Infusion Cannula. Straight

