

Vitreoretinal Cannulas and Vitreoretinal Handles

Silicone Brush Cannulas

Code	Gauge and size	Description	Tip
M4449	25G x 1 $\frac{3}{8}$ " 0.5 x 35mm	Silicone Brush. Silicone extends 0.7mm beyond end of cannula	
M4459	23G x 1 $\frac{1}{4}$ " 0.64 x 32mm	Silicone Brush. Silicone extends 0.7mm beyond end of cannula	
M4475	20G x 1 $\frac{1}{4}$ " 0.9 x 32mm	Silicone Brush. Silicone extends 1.25mm beyond end of cannula	
M4475A	20G x 1 $\frac{1}{4}$ " 0.9 x 32mm	Silicone Brush. Silicone extends 0.6mm beyond end of cannula	

Soft Tip Cannulas

Code	Gauge and size	Description	Tip
M4456	25G x 1 $\frac{1}{4}$ " 0.5 x 32mm	Soft Tip. Silicone extends 1mm beyond end of cannula	
M4458	23G x 1 $\frac{1}{4}$ " 0.64 x 32mm	Soft Tip. Silicone extends 2mm beyond end of cannula	
M4457	20G x 1 $\frac{1}{4}$ " 0.9 x 32mm	Soft Tip. Silicone extends 1mm beyond end of cannula	

Backflush Handles and Combination Sets

Code	Description
M4489A	'Active/Passive' Backflush Handle with finger control. Allows conversion from active to passive aspiration.
M4504 - set	25G x 1 $\frac{1}{4}$ " 0.5 x 32mm 'Active/Passive' Backflush Handle with 25G Soft Tip Cannula
M4530 - set	23G x 1 $\frac{1}{4}$ " 0.64 x 32mm 'Active/Passive' Backflush Handle with 23G Soft Tip Cannula
M4531 - set	20G x 1 $\frac{1}{4}$ " 0.9 x 32mm 'Active/Passive' Backflush Handle with 20G Soft Tip Cannula

Vitreoretinal Heavy Liquid Infusion Handles

Code	Description
M4498	Heavy Liquid Infusion Handle with finger control and a luer lok connection

Vitreoretinal Handles

Backflush Handles & Combination Sets

Code	Gauge and size	Description	Tip
M4489		Backflush Handle for Passive Aspiration with finger control and a luer lok cannula connection	
M4501 - set	20G x 1¼"	0.9 x 32mm Backflush Handle for Passive Aspiration with 20G Soft Tip Cannula set	Soft Tip
M4502 - set	23G x 1¼"	0.64 x 32mm Backflush Handle for Passive Aspiration with 23G Soft Tip Cannula set	
M4503 - set	25G x 1¼"	0.5 x 32mm Backflush Handle for Passive Aspiration with 25G Soft Tip Cannula set	
M4491 - set	20G x 1¼"	0.9 x 32mm Backflush Handle for Passive Aspiration with 20G Silicone Brush set	Silicone Brush
M4492 - set	23G x 1¼"	0.64 x 32mm Backflush Handle for Passive Aspiration with 23G Silicone Brush set	
M4496 - set	20G x 1¼"	0.9 x 32mm Backflush Handle for Passive Aspiration with 20G Aspiration Cannula set	Aspiration Cannula
M4497 - set	23G x 1¼"	0.64 x 32mm Backflush Handle for Passive Aspiration with 23G Aspiration Cannula set	
M4500 - set	25G x 1¼"	0.5 x 32mm Backflush Handle for Passive Aspiration with 25G Aspiration Cannula set	

Charles Flute Handles and Combination Sets

Code	Gauge and size	Description	Tip
M4493		'Charles' Flute Handle	
M4494 - set	20G x 1¼"	0.9 x 32mm 'Charles' Flute Handle with 20G Aspiration Cannula set	Aspiration Cannula
M4495 - set	23G x 1¼"	0.64 x 32mm 'Charles' Flute Handle with 23G Aspiration Cannula set	