

NEW

STERIMEDIX[®]

Silkann[®]
aesthetics

*aesthetic
cannula
innovations*




NEW

Silkann[®]
aesthetics

aesthetic cannula innovations

from Sterimedix - a world leader in single-use cannulas

*Introducing a new range of
flexible aesthetic cannulas -
designed with your patients'
comfort in mind*



Photograph of cannula tip
at 20x magnification

Sterimedix is a manufacturer of single use surgical cannula products based in the UK. With over 20 year's experience of manufacturing devices, primarily for ophthalmology, the Company has a rare depth of expertise in new product development.

The requirements of eye surgery place a great responsibility on the manufacturer to maintain the highest standards of quality. The regulatory environment in which the Company operates demands the best manufacturing and control standards.

Silkann is a new and complete range of aesthetic cannulas designed and manufactured to the highest levels of quality and consistency using state of the art manufacturing techniques.

- A new and unique screw type hub manufactured using polycarbonate to reduce the risk of cannula detachment.
- Indicator on the hub to show orientation of the port during injection.
- The best port quality and consistency of port position, to remove risks during procedures.
- Ports designed to optimise the flow and performance of fillers whilst reducing injection forces.
- The highest quality of packaging to maintain the integrity of the cannulas at all times.

For more information call: +44 (0)1527 405860 or visit: www.silkann

Code	Size	Description	Qty per box
------	------	-------------	-------------

18 GAUGE




AE 1701	18g UTW x 70mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40
---------	----------------	--	----

22 GAUGE




AE 1731	22g TW x 70mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40
AE 1732	22g TW x 60mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40
AE 1733	22g TW x 50mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40

23 GAUGE




AE 1740	23g UTW x 70mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40
AE 1742	23g UTW x 50mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40

25 GAUGE



AE 1750	25g UTW x 50mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40
AE 1751	25g UTW x 40mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40

27 GAUGE



AE 1761	27g UTW x 40mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40
---------	----------------	--	----



AE 1861	27g UTW x 40mm + 27g x 12.7mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port, packed with Sharp Needle	40
---------	----------------------------------	---	----

30 GAUGE



AE 1770	30g UTW x 30mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port	40
---------	----------------	--	----



AE 1870	30g UTW x 30mm + 30g x 12.7mm	Straight Dermal Cannula with Screw Type Hub, Closed End & Single Side Port, packed with Sharp Needle	40
---------	----------------------------------	---	----

Silkann[®] aesthetics

aesthetic cannula innovations



Sterimedix is based in Redditch, a large town just south of Birmingham, in the heart of the traditional needle making area of the UK. The first single use cannulas for ophthalmology were produced in Redditch over 40 years ago and the related skills remain in the area to this day. With a programme of investment in new technologies and automation, Sterimedix has been able to strengthen its position as a leader in the design, manufacture and sale of precision single use surgical devices.



With approvals to sell its products in all major markets world-wide, Sterimedix has developed strong relationships with both clinicians and other companies working in the field giving it a unique insight to their needs. Following a long history of manufacturing delivery devices for viscous products for ophthalmology, Sterimedix has now launched a new and unique range of cannulas designed specifically for use with aesthetic fillers.



Sterimedix Limited

1 Madeley Road, North Moons Moat, Redditch, Worcestershire B98 9NB UK

T : +44 (0)1527 405860 F : +44 (0)1527 501491

e : sales@sterimedix.com

www.silkann.com

STERIMEDIX[®]

is a registered trademark of Sterimedix Limited